



# NAUTILUS

SURGICAL STAPLERS



At Surgicaleuro Medical Devices Distributor, Ltd, patient healthcare is at the core of our everyday work.

We fully understand the challenging financial environment that healthcare gives face every day. So we work to provide a complete portfolio with the latest surgical technology and top quality products at an affordable cost thus helping to improve the life of more patients.

*“First Be Best Then Be First” – Grant Tinker*

Although we are a young company, we have at our heart more than 4 Years of expertise in the healthcare industry and in particular in the surgical field. All our staff have built a strong and solid reputation through the years and have a track record of excellence, business ethics, and commitment.

Only when we are the best can we reach first.

Our partners are a part of the company.

We have built a strong network of partners from suppliers to distributors that we consider as being a part of our company and build together the best solutions to be delivered and provided to patients worldwide.

Surgicaleuro Medical Devices Distributor, Ltd  
Your health. Our commitment.



# ABOUT



www.imq.it

CERTIFICATO N.  
CERTIFICATE N. 0517.2019

SI CERTIFICA CHE IL SISTEMA QUALITÀ DI  
WE HEREBY CERTIFY THAT THE QUALITY SYSTEM OPERATED BY  
**SURGICALEURO MEDICAL DEVICES  
DISTRIBUTOR, UNIP. LDA**

RUA CAROLINA MICHAELIS VASCONCELOS 30 CV, D1ª - 1500-145 LISBOA PORTUGAL

UNITÀ OPERATIVE / OPERATIVE UNITS

Vedere gli Allegati per le Unità Operative (n° 2 allegati)  
View the Annexes for the Operative Units (n° 2 annexes)

E' CONFORME ALLA NORMA / IS IN COMPLIANCE WITH THE STANDARD

**ISO 13485:2016**

PER LE SEGUENTI ATTIVITÀ / FOR THE FOLLOWING ACTIVITIES

Management of production and sale of disposable surgical  
staplers and accessories. Distribution of non-active and  
non-invasive medical devices for breathing and disposable gloves

Ulteriori informazioni riguardanti l'applicabilità dei requisiti ISO 13485:2016 possono essere ottenute consultando l'organizzazione  
Further clarifications regarding the applicability of ISO 13485:2016 requirements may be obtained by consulting the organization

IL PRESENTE CERTIFICATO E' SOGGETTO AL RISPETTO DEL  
REGOLAMENTO PER LA CERTIFICAZIONE DEI SISTEMI DI GESTIONE  
THE USE AND THE VALIDITY OF THE CERTIFICATE SHALL SATISFY THE  
REQUIREMENTS OF THE RULES FOR CERTIFICATION OF MANAGEMENT SYSTEMS



IONet, the association of the world's first class  
certification bodies, is the largest provider of Management  
System Certification in the world.  
IONet is composed of more than 20 bodies and counts  
over 100 subsidiaries all over the globe.



## EC CERTIFICATE

Certificate No 2026/MDD  
**Full Quality Assurance System Approval Certificate**

On the basis of our examination carried out according to Annex II, excluding section 4, of the Directive  
93/42/EEC and its revised version, we hereby certify that:

### Surgicaleuro Medical Devices Distributor Unip. Lda

1500-145 Lisboa - Rua Carolina Michaelis Vasconcelos 30 CV, D1ª (PRT) - Portugal  
manages in the factory of:

2700-391 Amadora - Parque Meramar VI, Armazén 1083 (PRT) - Portugal  
1500-145 Lisboa - Rua Carolina Michaelis Vasconcelos 30 CV, D1ª (PRT) - Portugal

a quality assurance system ensuring the conformity of the following products:

**Disposable surgical staplers and accessories**  
series and type refs in the Annex

with the relevant essential requirements of the aforementioned directive (from design to final inspection  
and testing) and it is subject to surveillance as specified in section 5 of Annex II. For class III devices, this  
certificate is valid only with the relevant EC Design-Examination Certificate of Annex II.4.

Reference to IMQ files Nos:  
DM18-0032937-01.

**This Approval Certificate is issued by IMQ S.p.A. as Notified Body for the Directive 93/42/EEC and its  
revised version. Notified Body notified to European Commission under number: 0051.**


Date:


2019-03-13

IMQ

# INDEX

Endoscopic Stapler	4
Circular Stapler	7
PPH Stapler	9
Linear Cutter	12
Linear Stapler	14





• Four kinds of cartridge:  
White, Blue, Gold and Green.  
Suitable for different  
thickness of tissue.

### Disposable Endoscopic Cutter Stapler and Cartridge


Model Stapler	Basic Size LV	Model Cartridge	Effetive Lenght	Staple Height		Staple Qty	Trocar
	Basic Size +-10		Basic Size +- 2	Code (color)	Height +- 2		
SE16-02-001	516	SE16-02-005	48.5	M Grey	2.0	70	12
SE16-02-002	616	SE16-02-006		W White	2.5		
		SE16-02-007		B Blue	3.5		
		SE16-02-008		D Gold	3.8		
		SE16-02-009		G Green	4.1		
SE16-02-003	528	SE16-02-010	60	W White	2.5	88	
SE16-02-004	628	SE16-02-011		B Blue	3.5		
		SE16-02-012		D Gold	3.8		
		SE16-02-013		G Green	4.1		

# Endoscopic Stapler 45 Premium

- 1 Adaptive joint head can swing within 90 degrees. Adjustment method is flexible, it can be operated by single hand.
- 2 Manually reset the button, the operation is safer.
- 3 Indicate the direction of the blade walking ensure safe operation without accident.
- 4 3+1 firing model, manual fi- ring back after 3 step firing, the stapler is more stable and safer.

Applicable to the tissue transection, resection and anastomosis in thoracic surgery, abdominal surgery and gynecologic endoscopy.





• Four kinds of cartridge:  
White, Blue, Gold and Green.  
Suitable for different  
thickness of tissue.

### Disposable Endoscopic Cutter Stapler and Cartridge

Model Stapler	Basic Size LY		Model Cartridge	Effetive Lenght L		Staple Height			Staple Qty	Trocar
	Basic Size	Limit Deviation		Basic Size	Limit Deviation	Code (color)	Height +- 2	Limit Deviation		
SE16-03-001	532.5	+- 10	SE16-02-006	48.5	+-2	W White	2.5	+- 0.2	70	12
			SE16-02-007			B Blue	3.5			
			SE16-02-008			D Gold	3.8			
			SE16-02-009			G Green	4.1			
SE16-03-002	544.5	+- 10	SE16-02-010	60	+-2	W White	2.5	+- 0.2	88	12
			SE16-02-011			B Blue	3.5			
			SE16-02-012			D Gold	3.8			
			SE16-02-013			G Green	4.1			


## Endoscopic Stapler R

- 1 Unique slot design for the front, which is convenient for the target organization.
- 2 Manual return knob - when the blade has no automatic repositioning, it is required to ensure safety.
- 3 3 step firing indication.

Applicable to the tissue transection, resection and anastomosis in thoracic surgery, abdominal surgery and gynecologic endoscopy.

The staple is the same when firing and firing is more smooth.





Model Stapler	Type	Length (mm)	Allowance
SE16-01-002	Standard	377	+ - 10
SE16-01-001	Short	277	
SE16-01-003	Long	477	+ - 10

**Disposable Endoscopic Cutter Stapler and Cartridge RS31\***

Model Cartridge	Color	Suture Length +- 5	Cutting Length +-5	Staple Height (mm)	Staple Width
SE16-01-011	White	45	40	2.5	
SE16-01-012	Blue			3.5	
SE16-01-013	Dark			4.1	
SE16-01-014	Green			4.8	
SE16-01-015	White	60	55	2.5	
SE16-01-016	Blue			3.5	
SE16-01-017	Dark			4.1	
SE16-01-018	Green			4.8	

\* RFS31 - Reinforced Staple 3 line

# Endoscopic Stapler 45


- 1 The joint head can swing 88° from left to right and rotate 360°.
- 2 The rotary knob of the joint is rotated and smooth, steady and non sloshing.
- 3 A more convenient cartridge load-off button.

Used in the tissue transection, resection and anastomosis of endoscopic surgery for visiting surgeon, gynecology, pediatrics, chest surgery.

Better blood supply, better tissue thickness adaptability, less tissue crush injury.

Fixed anvil design: When the anvil is closed to the jaw, the active side is on the side of the cartridge, the anvil side stay still, easy to pass through the narrow space, it is suitable for careful surgery and can protect tissue and blood vessels.






## Circular 3

- 1 Full firing position.
- 2 B-type staple with slightly in-curved feet make the anastomotic a higher tensile strength.
- 3 Three row titanium staples around can effectively close the tissue. The sharp circular knife can effectively cut.
- 4 The three row design is more reliable and the outer stapler is 0,2mm lower than the inner staple to insure good blood supply.
- 5 Safety lock button

Dedicated to the digestive tract reconstructive surgery for end-to-end, end-to-side and side-to-side anastomosis.

### Three Rows Series of Disposable Circular Stapler

Model Stapler	Color	Suture Rows	Outer Diameter mm	Cutting Diameter mm	Staple Height (mm)	Staple Width (mm)	Staple Qty	Staple Closure Height (mm)
SE16-04-001	Yellow Cartridge	3	21.2	12.4	4.4	3.7	27	2
SE16-04-002			24.1	15.0		3.4	33	
SE16-04-003			25.1	16.4		3.8		
SE16-04-004			26.1	17.0	4.7	3.9	36	
SE16-04-005			29.6	20.3			42	
SE16-04-006			31.0	21.0	4.9	3.9		
SE16-04-007			32.1	22.0			42	



## Circular

1 Safety lock button.

Applied to throughout the alimentary tract for end-to-end, end-to-side anastomose.

Large cutting diameter.

Adjustable thickness of tissue, meets all kinds of need.

Quicker opening and closing mechanism make it more comfortable to use.


Knife bunting seat without lip, it is convenient for instruments exit safely after firing.

### Disposable Circular Stapler

Model Stapler	Staple	Outer Diameter	Cutting Diameter mm	Width (mm)	Height (mm)
SE16-05-001	16	19.6	11.8	4.0	4.5
SE16-05-002		21.6	12.4		
SE16-05-003	18	24.3	15.0		
SE16-05-004	20	25.5	16.4		4.8
SE16-05-005		26.0	17.0		
SE16-05-006	24	29.7	20.7		5.0
SE16-05-007	26	31.9	23.1		







# PPH 3 Hemo

- 1 Full firing position, make the surgery safety.
- 2 Patent of three row, the effect of hemostasis is better and reduced the bleeding problem.
- 3 Dedicated to anal hemorrhoidectomy and excision of rectal prolapse.
- 4 Safety lock button.

## Three Rows Series of Disposable Circular Stapler

Model Stapler	Color	Suture Rows	Outer Diameter mm	Cutting Diameter mm	Staple Height (mm)	Staple Width (mm)	Staple Qty	Staple Closure Height (mm)
SE16-06-001	Yellow Cartridge	3	33.5	23.1	3.4	3.8	48	0.75-1.6
SE16-06-002			35.0	24.8	4.1	3.9	48	0.75-1.6





**Series of Disposable Anoscope \***

- Suitable for the examination of the lesions of the rectum, anus and so on.

\*Sold separately.



**Disposable Circular Stapler**

Model Stapler	Staple	Outer Diameter	Width (mm)	Height (mm)
SE16-08-001	32	34	4.0	4.3

**Visible PPH**

Suitable for the resection of anal hemorrhoids and internalrectal prolapse.





Circular Anal Dialator



Purse String Anoscope



Suture Threader



### Disposable Circular Stapler

Model Stapler	Staple	Outer Diameter	Width (mm)	Height (mm)
SE16-07-001	32	32.5	4.0	3.8
SE16-07-002		33.5		4.2

## PPH

Suitable for anal hemorrhoidectomy, rectal, & prolapse.





# Linear Cutter V1

Applied to the dissection, resection and suture of the digestive tract tissues of the stomach and intestines.

Blade in the cartridge.

## Stapler

Model Stapler	Basic Size	Cutting Line Length (mm)	Closure Height (mm)
SE16-09-001	255	60	1.0
SE16-09-002	275	80	1.5
SE16-09-003			2.0

## Cartridge

Model Cartridge	Color	Staple Height (mm)	Staple Qty
SE16-09-004	Blue	3.8	64
SE16-09-005			84
SE16-09-006	Yellow	4.3	



# Linear Cutter V2

Applied to the dissection, resection and suture of the digestive tract tissues of the stomach and intestines.



## Stapler

Model Stapler	Basic Size	Closure Height (mm)	Cutting Line Length (mm)
SE16-10-001	275	1.0	60
SE16-10-002	305	1.5	80
SE16-10-003	305	2.0	
SE16-10-004	335		100

## Cartridge

Model Cartridge	Color	Staple Height (mm)
SE16-10-005	Blue	3.8
SE16-10-006		
SE16-10-007	Yellow	4.3
SE16-10-008		



# Auto Linear 3R



Applicable to the digestive tract tissue residue end or incision closure.

The unique three row design, Unparalleled hemostatic effect, no need to reinforce the suture without special cases.

One-handed positioning, prevent spillover, easy operation and safer.

## Stapler

## Cartridge

Model Stapler	Colour	Suture Row	Lenght (mm)	Model Cartridge	Staple Height (mm)	Staple Width (mm)	Staple qty	Closure Height (mm)
SE16-13-001	Green 30	3	395	SE16-13-005	4.5	27	17	1.5-2.0
SE16-13-002	Green 45			SE16-13-006		26		
SE16-13-003	Green 60			SE16-13-007		35		
SE16-13-004	Green 90			SE16-13-008		53		





# Auto Linear

1 Tissue thickness  
1.5 - 2.5 mm

The closure of the stump or incision that is mainly applied to the digestive tract.

## Stapler

Model Stapler	Cutting Length (mm)	Staple Height (mm)	Length (mm)	Staple Qty
SE16-14-001	30	4.5	350	11
SE16-14-002	45			15
SE16-14-003	60			21
SE16-14-004	90			31





1

# Linear

1 Tissue thickness  
1.5 - 2.5 mm

Applied throughout the alimentary track for transection and resection of internal tissues.  
Four specifications meet different anastomosis lengths.

Closure height of staple can be adjusted according to the thickness of the tissue.

One of the spare cartridge can be used in the same operation.

## Linear

Model Stapler	Cutting Length (mm)	Staple Height (mm)	Length (mm)	Staple Qty
SE16-15-001	31.3	4.5	302	11
SE16-15-002	42.5			15
SE16-15-003	53.3			19
SE16-15-004	92.0			33







**SURGICALEURO**<sup>TM</sup>  
European medical devices



### **Headquarters:**

Rua Carolina Michaelis de Vasconcellos, 30 cv dta,  
1500-146 Lisbon, Portugal

### **Warehouse:**

Parque Meramar VI, Armazém 1083,  
2700-391 Amadora, Portugal

### **Swiss Office:**

Route de Barra 1c  
1974 Arbaz, Switzerland

**PT-Tel: 00351934731006**  
**CH-Tel.: 0041762119000**

**info@surgicaleuro.com**

**www.surgicaleuro.com**